

Activity Determination of Co-60 and Cs-137 Radiation Sources

Daniel Hackl

Executed at the TRIGA Center Atominstitut



Supervised by:

Dieter Hainz, BSc

Ing. Monika Veit

Dipl.-Ing. Dr. techn. Andreas Musilek

Motivation


- Inventory of historical radiation sources at the Triga Center Atominstitut
- Historically missing certificates
- Precise characterization required for proper disposal by NES GmbH
- Measurement methods/detectors used:
 - Dose rate meter (Thermo FH 40 G-10)
 - Geiger-Mueller-Counter
 - Activimeter (ISOMED 2010)
 - HPGe-Gammaspectrometer

Dose Rate Meter

- Measurement of the dose rate D at various measuring points
- Activity A can be determined using the square law of distance [1]:


$$A = \frac{(D_{Source} - D_{Bkgd}) \cdot \Gamma}{r^2} = \frac{D_{net} \cdot \Gamma}{r^2}$$

- Correction of the dose rate for Co-60 necessary
- Python script determines curve fit using the ODR method and returns activity A



Activimeter

- ISOMED 2010 from the manufacturer NuviaTech Healthcare for medical applications
- Multiple measurements in different measurement angles



Geiger-Mueller-Counter

- Measurement at a fixed distance of 8.5 cm
- Determination of detector efficiency using two reference sources:

| Isotope | Energy (keV) | I_γ | η (8.5 cm) |
|---------|--------------|------------|----------------------|
| Co-60 | 1252,86 | 0,9992 | $1,01 \cdot 10^{-4}$ |
| Cs-137 | 661,655 | 0,8505 | $2,09 \cdot 10^{-5}$ |

- Average activity (three measurements each) [2]:

$$\bar{A} = \frac{\bar{n}}{\eta(E_\gamma) \cdot t_{meas} \cdot I_\gamma} \quad \sigma_{\bar{A}} \approx \frac{\sqrt{\bar{n}}}{\sqrt{N} \cdot \eta(E_\gamma) \cdot t_{meas} \cdot I_\gamma}$$




Gammaspectroscopy

- High purity Germanium detector
- Co-60 has two γ -peaks, activity calculated with equation [3]:

$$A = \sum_i \frac{I_i(E_\gamma)}{t_{real} \cdot \eta(E_\gamma) \cdot R_t \cdot I_\gamma}$$

- Efficiency for different geometries and isotopes calculated with reference sources



Results

| Source | Calc. A | Dose Rate | Activ. | γ -Spectr. | Geiger-C. |
|-----------------|---------|-----------|--------|--------------------|----------------|
| Nr. 2 (Co-60) | 3,621 | 3,224 | 3,685 | - | 3,221 |
| Nr. 3 (Co-60) | 1,254 | 1,239 | 1,337 | - ¹ | - ² |
| Nr. 4 (Co-60) | 0,276 | 0,782 | 0,941 | 0,837 | 0,860 |
| Nr. 5 (Co-60) | 0,001 | - | - | 0,002 | - |
| Nr. 6 (Co-60) | 0,041 | - | 0,008 | 0,015 | - |
| T 15/4 (Cs-137) | - | - | 0,052 | 0,041 ³ | - |
| OTL 99 R004 | - | 0,544 | 0,518 | 0,508 | - |
| OTL 99 R005 | - | 0,534 | 0,535 | - ⁴ | - |
| FL 90/1975 | 0,110 | 0,097 | 0,100 | 0,096 | - |
| automess 6706 | 0,126 | 0,134 | 0,106 | 0,098 | - |
| T25-20 | - | 4,757 | 4,502 | - ⁵ | 4,584 |
| T 3/4 | 0,091 | - | 0,087 | 0,057 | - |
| Nr. 7 | 0,220 | 0,137 | 0,145 | 0,135 | - |
| Nr. 8 | 58,127 | 48,317 | 45,418 | 37,702 | - ⁶ |
| Nr. 30 | 0,316 | 0,291 | 0,309 | 0,262 | - |

Discussion

- **Dose Rate Meter:** Precise values for a wide range of activities
- **Activimeter:** Not consistent, high error if below minimum activity value
- **Geiger-Mueller-Counter:** Suitable for activities $> 0.8 \text{ MBq}$
- **γ -Spectroscopy:** Not consistent (efficiency), not suitable for high activities, best results for activity range $1 \text{ kBq} - 550 \text{ kBq}$

References

- [1] M. Tschurlovits, A. Leitner und G. Daverda, »Dose Rate Constants for New Dose Quantities,« Radiation Protection Dosimetry, Jg. 42, Nr. 2, p. 77-82, Aug. 1992.
- [2] N. Tsoulfanidis und S. Landsberger, Measurement and Detection of Radiation. CRC Press, Taylor & Francis Group, 2015, ISBN: 9781482215496.
- [3] »Praktikum aus Neutronenphysik: Aktivierung von Gold- und Indiumfolien,« 2022.